

WOUND MANAGEMENT FORMULARY 2011

PRODUCT INFORMATION

**FOR CARE HOMES AND EXTERNAL
PARTNERS**

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Product Generic Grouping	1st Line Products	2nd Line Products
Hydrogel Hydrogel Sheet	ActivHeal Hydrogel (Advanced Medical Solutions) Actiform Cool hydrogel sheet (Activa).	Askina Gel (B-Braun) Nil
Hydrocolloid	Duoderm Extra thin (Convatec)	Nil
Foams for Low Exudate (Adhesive & Non-adhesive)	PolyMem (Aspen Medical) (Not shaped products)	Nil
Foams for Moderate to High Exudate (Adhesive & Non-adhesive)	Nil	Tielle Range (Systagenix) (Not shaped products)
Hydrofibre	ActivHeal Aquafiber (Advanced Medical Solutions)	Aquacel (Convatec)
Foams – shapes (heels)	Tegaderm Heel (3M)	Nil
Silver	Silvercel Non-Adherent (Systagenix) Actisorb Silver 220 (Systagenix)	Nil Nil
Honey	Actilite, Activon Tube, Activon Tulle, Algivon, (Advancis Medical)	Nil
Bio-enzyme	Flaminal Hydro & Flaminal Forte (Ark Therapeutics)	Nil

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Increased Absorption	Eclipse (Advancis Medical)	Exu-dry (for specialist/palliative care use only)
Wound Contact Layer	NA Ultra (Systagenix)	Adaptic Touch (Systagenix)
Silicone Foam (Border and Non-Bordered)	Mepilex (Molnycke)	Allevyn Gentle
Charcoal- odour absorbing dressings.	Clinisorb (Clinimed)	Carboflex (Convatec) for specialist/palliative care use only
	Special Access	
	1st Line Products	Advice
Iodine	Inadine (Systagenix) Iodoflex, (S&N)	Do not use in thyroid disease and renal impairment. Use for stubborn infections especially MRSA.
Larvae	Bio-Monde (formerly Zoobiotic) via GP/Pharmacy	
	Via Specific Specialties	
<u>Dermatology</u> Paraffin Gauze Antibiotic Gel	Jelonet (S&N) Metronidazole gel (via pharmacy) Flamazine (S&N) (via pharmacy)	Plastics use also
<u>Podiatry</u> Silver Foam Foam Boot	Biatain Foam Biatain Ag Foam (Coloplast) Askina Heel Foam Kerraboot (Ark Therapeutics)	

Disclaimer: For more in depth product information please refer to company product leaflet.

**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

Tissue Viability Topical Negative Pressure (TNP)	Venturi (Talley)	
Protease Modulator	Start (Urgo) Promogram Prisma	
Sundries		
Sticky tapes Adhesives Cleanser	Clinipore Mefix Prontosan Solution	
Bandages Retention products	To be confirmed @ Nov 2011 Comifast Range	
Barrier creams etc	Cavilon Barrier Cream (3M) Cavilon barrier Film (3M) Zinc paste or equivalent	

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category	Hydrogel Sheet
Product	Actiform Cool
Company	Activa Healthcare
Indications for Use	<p>Hydrates dry tissue (sloughy or necrotic). Provides pain relief through its cooling effect. May be used-</p> <ul style="list-style-type: none"> • at all stages of wound healing • May be used with compression bandaging for moderately to heavily exuding wounds. • May be used on contaminated and colonised wounds but systemic antibiotics would be required. • Haemostasis should be achieved prior to use on granulating wounds with haemoserous (bloody) fluid. • Provides relief for irritated, sore skin post-radiotherapy treatment.
Contra-Indications for Use	<p>Not to be used as a covering on deep, narrow cavities or sinuses. The wound should be checked regularly as the absorption process may be quite rapid and the wound may become drier than expected. Information on the product's "liner" should be read before use. Contains phenoxyethylene as a preservative no noted sensitivities</p>
Application Guidance	<p>The sheets can be cut to size or placed side by side for larger areas. The sheets can be removed in one piece if the gel has not yet formed. If the gel has formed, it can be scooped up or washed/soaked off gently with water. The dressing should be changed when the gel becomes cloudy or opaque or at the first signs of strike through.</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	084159 (304-8352)
10cm x 15cm	093205 (304-8360)
5cm x 6.5cm	117123 (315-5553)

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Hydrogel
Product	ACTIVHEAL HYDROGEL Contains 85% water and a collection of polymer chains that are water insoluble
Company	Advanced Medical Solutions
Indications for Use	Mainly used to donate liquid to - <ul style="list-style-type: none"> • Dry slough or necrotic wound • Granulating wounds • Nil to Low exuding wounds
Average Wear Time	Dependant on exudate and leakage. Can be left in place up to 5 days. Monitor more closely in PVD and Diabetes
Contra-Indications for Use	<ul style="list-style-type: none"> • Do not use on infected wounds • Do not use on patients with known sensitivity to the gel or one of its components • Do not use prior to Larvae therapy
Application Guidance	<ul style="list-style-type: none"> • Apply directly into the wound at a minimum layer of 5mm • Cover with NA Ultra initially • Then cover with a secondary dressing. • Remove by irrigation with normal saline. • Single use only - discard remnants
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
15 gm	SKU 130474 /1007419

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 2nd Line Choice	Hydrogel
Product	ASKINA GEL Based on unique formulation which provides more hydration to the wound and demonstrates an outstanding stability.
Company	B-Braun
Indications for Use	<ul style="list-style-type: none"> • Dry slough or necrotic wound • Nil to Low exuding wounds • Granulating wounds • Narrow wounds and sinuses
Average Wear Time	Dependant on exudate and leakage. Can be left in place up to 5 days. Monitor more closely in PVD and Diabetes
Contra-Indications for Use	<ul style="list-style-type: none"> • Do not use on infected wounds • Do not use on patients with known sensitivity to the gel or one of its components • Do not use on third degree burns.
Application Guidance	Product tip/cannula allows easy access to deeper wounds such as tracks and sinuses. <ul style="list-style-type: none"> • Apply directly into the wound at a minimum layer of 5mm • Cover with NA Ultra initially • Then cover with a secondary dressing. • Remove by irrigation with normal saline. • Single use only - discard remnants
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
15gm	134588

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st line Choice	Hydrocolloid
Product	Duoderm Extra Thin Hydrocolloid dressing with waterproof outer layer which is impermeable to exudates and micro-organisms The hydrocolloid component of the dressing interacts with the wound exudates to produce a cohesive gel which enables removal with little or no damage to newly formed tissue.
Company	Convatec
Indications for Use	Minimal or lightly exudating wounds and minor skin injuries such as superficial pressure sores and sutured wounds. <ul style="list-style-type: none"> • To rehydrate necrosis and dry slough • Prophylactic for prevention of friction damage due to fragile skin, high risk patients. • Able to bathe or shower with dressing in position.
Average wear-time	Dependant on exudate and leakage. Can be left in place up to 7 days. Monitor more closely in PVD and Diabetes.
Contra-Indications for Use	Not ideal for infected wounds. <ul style="list-style-type: none"> • Use of Duoderm dressings in the presence of anaerobic infection is not recommended. • Colonisation of chronic wound is common and not contraindicated to the use of Duoderm dressings. • Where systemic or local infection develops during the use of the dressing, appropriate adjunctive therapy should be initiated.
Application Guidance	When applying the dressing ensure minimum 2 cm overlap from wound margin to ensure good adhesion to surrounding skin.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	031658
15cm x 15cm	031665
5cm x 10cm	034857
5cm x 20cm	S164
7.5cm x 7.5cm	034840
9cm x 15cm	S171
9cm x 25cm Jubilee/Ortho use	078677
9cm x 35cm Jubilee/Ortho use	078684

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	FOAMS for LOW EXUDATE
Product	POLYMEM PolyMem dressings are thin polyurethane foam membrane dressings containing a non-toxic cleanser (surfactant) and a moisturiser (glycerol), with a semi-permeable film backing. The unique properties of a PolyMem dressing also help to reduce inflammation, promoting comfort and pain reduction while the dressing is worn.
Company	Aspen Medical
Indications for Use	<ul style="list-style-type: none"> • All wounds • Fragile skin • Low to moderate exuding wounds • Minimises procedural pain and trauma • May increase exudates initially
Average Wear Time	Depends on exudate levels
Contra-Indications for Use	<ul style="list-style-type: none"> • Discontinue use if any sensitivities occur. • Do not use if sensitive to any glycerol
Application Guidance	<ul style="list-style-type: none"> • Primary dressing only • No cleansing required unless debris/contaminants visible • Use barrier film to surrounding tissue • If using non-adhesive - do not cover with film dressing.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
Non-Adhesives	
8cm x 8cm	5033
13cm x 13cm	5055
17cm x 19cm	5077
Roll 20cm x 60cm	5824
PolyMem WIC 8 x 8 cm for shallow cavities	5733
Adhesives	203
5 x 5cm dot (2.5 x 2.5cm pad)	
Oval 5 x 7.6cm (2.5 x 5cm pad)	8023
Oval 8.8 x 12.7cm (5 x 7.6cm pad)	8053
15 x 15cm island (9 x 9cm pad)	606
Oval 16.5 x 20.9cm (10.1 x 14.6cm pad)	8086

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Moderate to High Exuding Wounds Foam
Product	TIELLE (Adherent) TIELLE PLUS (Adherent) TIELLE EXTRA (Non-adherent)
Company	Systagenix
Indications for Use	<ul style="list-style-type: none"> Assess exudate levels and skin condition and choose appropriate product Moderate to high exuding wounds All types of wounds
Average Wear Time	<ul style="list-style-type: none"> Depends on exudate levels
Contra-Indications for Use	<ul style="list-style-type: none"> Third degree burns Active vasculitis
Application Guidance	Use as a primary or secondary dressing. If using non-adherent secure with appropriate retention dressings or bandage Use tap water along the product edges for atraumatic removal
PCOD Ordering information	Product brand name, size and quantity required

Range	Available Sizes	Procurement Codes
Tielle Adhesive	11 x 11 cm	SKU 065981
	15 x 15 cm	SKU 135127
	15 x 20 cm	SKU 065998
Tielle Plus Adhesive	11 x 11 cm	SKU 061952
	15 x 15 cm	SKU 134106
	15 x 20 cm	SKU 135134
Tielle Extra Non-adherent	11 x 11 cm	MTP301
	15 x 15 cm	MTP305
	15 x 20 cm	MTP302

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Hydrofibre
Product	AQUAHEAL AQUAFIBER
Company	Advanced Medical Solutions
Indications for Use	<ul style="list-style-type: none"> • Moderate to heavily exuding wounds • All wound types • To control minor bleeding in superficial wounds
Average Wear Time	<ul style="list-style-type: none"> • Remove when clinically indicated i.e. leakage, excessive bleeding, or suspicion of infection. • Can be left in place for up to 7 days
Contra-Indications for Use	<ul style="list-style-type: none"> • Sensitivity to this dressing or one of its components. • Do NOT use in sinuses and tracks
Application Guidance	Choose size slightly larger than wound. Apply direct to wound, loosely to fill deeper wounds. Apply secondary dressing
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
5x5cm	10009146
10x10cm	10009145
15x15cm	10009147

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 2nd Line Choice	Hydrofibre
Product	AQUACEL HYDROFIBRE
Company	Convatec
Indications for Use	<p>Aquacel Hydrofibre dressing may be used for the management of exuding wounds:-</p> <ul style="list-style-type: none"> • Leg ulcers, pressure ulcers and diabetic ulcers. • Surgical wounds • Partial thickness burns • Traumatic wounds • Exudate absorption in oncology wounds • Use in conjunction with systemic antibiotic therapy if infection noted.
Average Wear Time	<ul style="list-style-type: none"> • Remove when clinically indicated i.e. leakage, excessive bleeding, or suspicion of infection. • Can be left in place for up to 7 days.
Contra-Indications for Use	<ul style="list-style-type: none"> • Sensitivity to this dressing or one of its components. • Do NOT use in sinuses and tracks
Application Guidance	Dressing should overlap 1 cm on to the skin surrounding the wound. Only pack deep wounds up to 80%, as Aquacel will expand. Apply secondary dressing
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
5cm x 5 cm	031726
10cm x 10cm	031733
15cm x 15cm	031740

Category	Shaped Foams
-----------------	---------------------

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

1st Line Choice	Polyurethane Foam Dressing
Product	TEGADERM FOAM HEEL
Company	3M
Indications for Use	<ul style="list-style-type: none"> • Moderate to high exuding wounds. • Provide fast absorption of exudates, reducing the risk of maceration and the potential for skin breakdown
Average Wear Time	<ul style="list-style-type: none"> • Dependant on wound exudate, change when exudate is noted near the outer part of the grid. • Can be in place up to 7 days
Contra-Indications for Use	<ul style="list-style-type: none"> • Not for use in dry wounds. • Do not use on dry, scaly skin or excoriated skin. • Discontinue use if redness or sensitivities noted
Application Guidance	<p>3M Foam Dressing (non-adhesive) is designed with a white, non – adherent wound contact side and an outer, barrier film side with a printed blue grid and 3M logo (apply grid facing up)</p> <p>Adhesive – remove backing and smooth gently over wound and surrounding skin</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
Heel 14cm x 14cm	034093

Category	Honey Oil
-----------------	------------------

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

1st Line Choice	Low Adherence
Product	ACTILITE Knitted viscose for low adherence and coated with anti-bacterial Activon+ which contains Manuka honey oil
Company	Advancis Medical
Indications for Use	All acute and chronic wounds at the epithelialisation stage, <ul style="list-style-type: none"> • Reduces odour • Antibacterial • Use as substitute for Inadine as more patient friendly with minimal risk of complications
Average Wear Time	Dependent on exudate level
Contra-Indications for Use	Known allergy to bee venom. Blood sugar levels should be monitored in patients with diabetes. Increased pain can occur due to osmotic pulling if unacceptable, remove dressing.
Application Guidance	The dressing is easy to apply, with honey being naturally self-adherent it will gently adhere to surrounding dry skin but will not adhere to a wound bed. Apply a secondary dressing of choice; this can be an absorbent dressing or semi-permeable film dressing depending on the nature of the wound.
PCOD Ordering information	Product brand name, size and quantity required

Sizes	Procurement Codes
10 x 10cm	138104
10 x 20cm	CR3852

Category	Antimicrobial
-----------------	----------------------

Disclaimer: For more in depth product information please refer to company product leaflet.

**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

1st Line Choice	Honey Ointment
Product	ACTIVON TUBE
Company	Advancis Medical
Indications for Use	All acute and chronic wounds, <ul style="list-style-type: none"> • Sloughy, necrotic and malodorous. • Reduces odour • Antibacterial
Average Wear Time	<ul style="list-style-type: none"> • Dependent on exudate level
Contra-Indications for Use	<ul style="list-style-type: none"> • Known allergy to bee venom. • Blood sugar levels should be monitored in patients with diabetes. • Increased pain can occur due to osmotic pulling if unacceptable, remove dressing.
Application Guidance	Apply tube directly to wound, cover with ActivonTulle or Algivon to maintain contact with wound bed. Observe for maceration of surrounding skin May also increase exudate levels, ensure use of adequate adsorbent secondary dressing.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
Activon Tube – 25g	113613

Category	Antimicrobial
-----------------	----------------------

Disclaimer: For more in depth product information please refer to company product leaflet.

**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

1st Line Choice	Honey Impregnated Gauze
Product	ACTIVON TULLE
Company	Advancis Medical
Indications for Use	All acute and chronic wounds, <ul style="list-style-type: none"> • All wound types – granular, sloughy, necrotic • Reduces odour • Antibacterial
Average Wear Time	<ul style="list-style-type: none"> • Dependent on exudate level
Contra-Indications for Use	<ul style="list-style-type: none"> • Known allergy to bee venom. • Blood sugar levels should be monitored in patients with diabetes. • Increased pain can occur due to osmotic pulling if unacceptable, remove dressing.
Application Guidance	Can be cut to size, Apply directly over wound. Layer can be unfolded and laid out for larger areas. Observe for maceration of surrounding skin May also increase exudate levels; ensure use of adequate adsorbent secondary dressing.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
Activon Tulle 10 x 10 cm	134052
Activon Tulle 5 x 5 cm	112494

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Antimicrobial Honey Impregnated Alginate
Product	ALGIVON
Company	Advancis Medical
Indications for Use	All acute and chronic wounds, <ul style="list-style-type: none"> • Use on dry necrosis to increase honey donation • Can also absorb exudate • Sloughy, necrotic and malodorous wounds. • Reduces odour • Antibacterial
Average Wear Time	<ul style="list-style-type: none"> • Dependent on exudate level
Contra-Indications for Use	<ul style="list-style-type: none"> • Known allergy to bee venom. • Blood sugar levels should be monitored in patients with diabetes. • Increased pain can occur due to osmotic pulling if unacceptable, remove dressing.
Application Guidance	Can be cut to size, Apply directly over wound. Observe for maceration of surrounding skin May also increase exudate levels; ensure use of adequate adsorbent secondary dressing.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	138098
5cm x 5cm	138081

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Antimicrobial Bio-Enzyme Hydrogel
Product	FLAMINAL FORTE
Company	Ark Therapeutics
Indications for Use	<ul style="list-style-type: none"> • Moderate to heavy exuding wounds. • Second degree burns, • Pressure ulcers • Cavities and sinuses
Average Wear Time	1 - 4 days depending on exudate levels.
Contra-Indications for Use	<ul style="list-style-type: none"> • Sensitivity reaction to alginate dressings • Sensitivity to polyethylene glycol (Macrogol)
Application Guidance	<p>Apply a thick layer of Flaminal about 5mm, avoiding the wound edge Apply NA Ultra over this to aid retention of gel then an appropriate secondary dressing. Leave in place whilst the gel structure remains intact Flaminal is self sterilising, so suitable for multi-use on a single patient</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
15g tubes	138067
50g tubes	PCOD use only

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Antimicrobial Bio-Enzyme Hydrogel
Product	FLAMINAL HYDRO
Company	Ark Therapeutics
Indications for Use	<ul style="list-style-type: none"> • Light exuding wounds. • Dermatosurgery, • Second degree burns • Pressure ulcers
Average Wear Time	<ul style="list-style-type: none"> • 1 - 4 days depending on exudates levels.
Contra-Indications for Use	<ul style="list-style-type: none"> • Sensitivity reaction to polyethylene glycol (Macrogol)
Application Guidance	<p>Apply a thick layer of Flaminal about 5mm, avoiding the wound edge Apply NA Ultra over this to aid retention of gel then an appropriate secondary dressing. Leave in place whilst the gel structure remains intact Flaminal is self sterilising, so suitable for multi-use on a single patient</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
15g tubes	138074
50g tubes	PCOD use only

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	High Absorbency
Product	ECLYPSE NON-ADHERENT
Company	Advancis Medical
Indications for Use	<ul style="list-style-type: none"> • High exuding wounds
Average Wear Time	<ul style="list-style-type: none"> • Wear time will depend on the level of exudate. • The dressing must be changed prior to reaching capacity.
Contra-Indications for Use	<ul style="list-style-type: none"> • Discontinue use if redness or sensitivities noted
Application Guidance	<p>Eclipse is placed face down on wound surface with beige backing uppermost.</p> <p>Due to the high absorbent capacity the dressing may need to be supported by an additional fixation e.g. a light support bandage or appropriate tape.</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	130757
15cm x 15cm	130740
20cm x 30cm	096268
60cm x 40 cm	138043
Boot 60 x 70 cm	CR4014

Disclaimer: For more in depth product information please refer to company product leaflet.

**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

Category Special Access	High Absorbent Palliative Use
Product	EXU-DRY
Company	Smith and Nephew
Indications for Use	<ul style="list-style-type: none"> • Very high exuding wounds
Average Wear Time	<ul style="list-style-type: none"> • Dependent on exudate levels
Contra-Indications for Use	<ul style="list-style-type: none"> • Known skin sensitivities to product
Application Guidance	<p>Select a size slightly larger than the wound. Place onto wound. The back side of the dressings are marked "Use other side against wound".</p> <p>Secure in place with appropriate tape around side edges (window pane) or Clinifast retention products or use a light bandage.</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
15cm x 23cm	5999006
23cm x 38cm	5999006
38cm x 46cm	5999018
38 cm x 61cm	5999024
Arm Dressing	5999LPA
Leg Dressing	5999LPL
Large Burn Vest	5999LV1

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Wound Contact Layer Lightly coated silicone gauze
Product	NA ULTRA
Company	Systagenix
Indications for Use	<ul style="list-style-type: none"> For use on all granulating wounds including burns, skin trauma, leg ulcers
Average Wear Time	<ul style="list-style-type: none"> Can be left up to 7 days
Contra-Indications for Use	<ul style="list-style-type: none"> Sensitivities to silicone.
Application Guidance	<p>Can be cut to size. Place over wound secure with secondary dressing depending on exudate levels Apply over gel and ointment based product to ensure good contact with wound.</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
19cm x 9.5cm	025367
9.5cm x 9.5cm	025374

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 2nd Line Choice	Wound Contact layer Silicone Coated Mesh
Product	ADAPTIC TOUCH
Company	Systagenix
Indications for Use	<ul style="list-style-type: none"> • Dry to highly exuding wounds • All wound types • Aids atraumatic removal • Reduces wound related pain on dressing changes
Average Wear Time	<ul style="list-style-type: none"> • Can be left in place up to 7 days
Contra-Indications for Use	<ul style="list-style-type: none"> • Known sensitivities to cellulose or silicone • Surgical implantation
Application Guidance	Choose appropriate size, can be cut to fit. Apply directly over the wound and cover with secondary dressing
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
5 x 7.6 cm	TCH501
7.6 x 11 cm	TCH502
12 x 15 cm	TCH503
20 x 32 cm	TCH504

Disclaimer: For more in depth product information please refer to company product leaflet.

**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

Category 1st Line Choice	Silicone Foam Adhesive & Non-Adhesive
Product	MEPILEX and MEPILEX BORDER
Company	Molnlycke Health Care
Indications for Use	<ul style="list-style-type: none"> • Fragile and sensitive skin • Pressure ulcers, • Leg and foot ulcers and • Traumatic wounds.
Average Wear Time	<ul style="list-style-type: none"> • Dependent on exudate levels
Contra-Indications for Use	<ul style="list-style-type: none"> • Discontinue use if reddening or sensitisation occurs. • Do NOT routinely use shaped products
Application Guidance	Prepare and clean the skin surrounding the wound area by removing excess moisture. Any excess hair should be clipped to ensure close approximation to the wound. Remove the protector material gradually and anchor the adhesive side of the dressing to the skin and ensure that the dressing is adhered all around the wound.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
Border 7cm x 7.5 cm	123636
Border 10cm x 12.5cm	123643
Border 10cm x 20cm	136698
Border 15cm x 17.5 cm	123650
Border 10cm x 30cm	117482
Non-Adhesive Range	
10cm x 11 cm	SKU 123599
11 x 20cm	294260
15 x 16cm	SKU 123612
20 x 21cm	SKU 123629

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 2nd Line Choice	Silicone Foam Adhesive & Non-Adhesive
Product	ALLEVYN GENTLE and ALLEVYN GENTLE BORDER
Company	Smith & Nephew
Indications for Use	<ul style="list-style-type: none"> • Chronic and acute exuding wounds, • shallow, granulating wounds. • full and partial thickness wounds
Average Wear Time	<ul style="list-style-type: none"> • Dressings can be left in place undisturbed for up to 7 days • Change when exudate is visible and approaches to within 1.5cm of the edge of the dressing
Contra-Indications for Use	Discontinue use if reddening or sensitisation occurs.
Application Guidance	<p>Prepare and clean the skin surrounding the wound area by removing excess moisture. Any excess hair should be clipped to ensure close approximation to the wound. Using aseptic technique remove the protector material and apply the adherent side to the wound ensuring good contact all the way around.</p> <p>For non bordered secure with stockinette retention product or bandage.</p>
PCOD Ordering information	Product brand name, size and quantity required

Range and Sizes Allevyn Gentle Border	Procurement Codes
7.5cm x 7.5cm	133819
10cm x 10cm	133826
12.5cm x 12.5cm	133833
17.5cm x 17.5cm	133840
Allevyn Gentle	
10cm x 10cm	66800248
10cm x 20cm	66800249
15cm x 15cm	66800250
20cm x 20cm	66800251
5cm x 5cm	66800247

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Charcoal Odour Absorbing
Product	CLINISORB
Company	Clinimed
Indications for Use	<ul style="list-style-type: none"> • For malodorous wounds, e.g. malignancy,
Average Wear Time	<ul style="list-style-type: none"> • Change when wet. • Can be left in place up to 7 days
Contra-Indications for Use	No known contraindications
Application Guidance	Can be used as both primary and secondary dressings. Place either side down to wound, May be cut to size.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	024728
10cm x 20cm	024735
15cm x 25cm	2315

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 2nd Line Choice	Charcoal Odour Absorbing
Product	CARBOFLEX Multi-layered odour absorbing dressing that will absorb exudate
Company	Convatec
Indications for Use	<ul style="list-style-type: none"> For exuding, malodorous wounds e.g. malignancy
Average Wear Time	<ul style="list-style-type: none"> Can be left up to 48 hours depending on exudate
Contra-Indications for Use	<ul style="list-style-type: none"> Known sensitivity to product components
Application Guidance	Select appropriate size dressing to overlap wound by 3 cm. Apply fibrous side of dressing to wound. Do not cut to size. Can be used as a secondary dressing over cavity wounds.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	117277
15cm x 20cm	034871
8cm x 15cm	049646

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	Iodine Gauze impregnated with iodine
Product	INADINE
Company	Systagenix
Indications for Use	<ul style="list-style-type: none"> • Minor skin trauma with stubborn infection • MRSA contaminated wounds
Average Wear Time	<ul style="list-style-type: none"> • Colour change indicates when dressing needs changing. • Medical supervision should be sought if using Inadine dressing for more than one week.
Contra-Indications for Use	<ul style="list-style-type: none"> • Do not use in known iodine hypersensitivity • Do not use in renal disease • Use with caution in thyroid disease • Consider use of Actilite as more patient – friendly.
Application Guidance	Remove one backing paper, peel off the dressing from the remaining backing paper and apply it directly to the wound. Cover and secure with a secondary dressing of choice
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
5cm x 5cm	039654
9.5cm x 9.5cm	061662

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	Iodine Impregnated
Product	IODOFLEX
Company	Smith & Nephew
Indications for Use	<ul style="list-style-type: none"> • Stubborn infection in chronic exuding wounds. • Iodoflex has dual action antimicrobial and desloughing properties.
Average Wear Time	<ul style="list-style-type: none"> • Dependant on the exudate levels. • Changes in colour indicate when dressing change is needed.
Contra-Indications for Use	<ul style="list-style-type: none"> • Do not use on dry necrotic tissue or on patients with a known sensitivity to any of its ingredients • Do not use on children, pregnant or lactating women or people with thyroid disorders or renal impairment
Application Guidance	<p>A single application of Iodoflex should not exceed 50g and not more than 150g in one week. The duration of treatment should not exceed 3 months and allow a one week break before restarting treatment to enable excretion of iodine.</p> <p>Remove carrier gauze on both sides of the dressing and apply to wound surface. The dressing can be moulded or cut to fit wound. Cover with an appropriate secondary dressing</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Via Pharmacy
5g	
10g	
17g	

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Larvae Sterile maggot
Product	BIOMONDE- LARVAE Unlicensed Product
Company	BioMonde (formerly Zoobiotics)
Indications for Use	<ul style="list-style-type: none"> • Debridement, • Infected rehydrated necrotic wounds , • Can potentially speed up healing times • Measure wounds to ensure correct size ordered.
Average Wear Time	Dressings need checked and turned daily. 1 dressing can last 5 days.
Contra-Indications for Use	Warfarin therapy Discontinue use if fresh bleeding noted.
Application Guidance	Larvae and sudocrem come boxed from company. Apply Sudocrem to surrounding skin Apply BioFoam pouch and cover with wound pad and light bandage or mefix to secure. Do not use films or foams, as larvae need oxygen to survive.
Acute Sites Ordering Information	Prescribe on medical Drug Kardex as once only -for each application Consultant to complete unlicensed medicine form Complete consent form
PCOD Ordering information	Product brand name, size and quantity required Add 'C' to end of below codes. Order before 12 midday for next day delivery. Can be delivered to pharmacy or health centre.

Sizes	Order Via Pharmacy Codes
2cm x 2cm	BFD 20
2.5cm x 4cm	BFD 25
5cm x 5cm	BFD 50
2.5cm x 15cm	BDF150
7cm x 7cm	BFD 70
7cm x 12cm	BFD 120
10cm x 10cm	BFD 100

Replace coded D with M if maintenance dose required

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	Paraffin Impregnated Gauze
Product	JELONET
Company	Smith & Nephew
Indications for Use	<ul style="list-style-type: none"> • Some Dermatology Conditions • Plastic Surgery
Average Wear Time	<ul style="list-style-type: none"> • Once – Twice daily
Contra-Indications for Use	<ul style="list-style-type: none"> • For specialist monitored use only. • Can cause trauma to newly formed granulation tissue.
Application Guidance	Remove backing and apply to wound using aseptic technique.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10cm x 10cm	017409
10cm x 40cm	027514
10cm x 7cm Tin	037438
15cm x 2cm Roll	070589
5cm x 5cm	070572

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category 1st Line Choice	Metronidazole Gel Excipients include hydroxybenzoates (parabens), propylene glycol.
Product	ANABACT GEL 0.75%
Company	
Indications for Use	<ul style="list-style-type: none"> • Malodorous fungating tumours. • Malodorous gravitational and decubitus ulcers.
Average Wear Time	1-2 times daily
Contra-Indications for Use	Skin irritation. Pregnancy. Allergy to – Metronidazole, Bronopol and hydroxybenzoic acid esters, hydroxyethylcellulose, propylene glycol,
Application Guidance	Apply the gel over the complete area and cover with a non adherent dressing.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
15 gm	via pharmacy
30 gm	via pharmacy

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	Silver Sulphadiazine 1% cream
Product	FLAMAZINE
Company	Smith and Nephew
Indications for Use	<ul style="list-style-type: none"> • prophylaxis and treatment of infection in • burn wounds • short-term treatment of infection in leg ulcers and pressure sores, • prophylaxis of infection in skin graft donor sites and extensive abrasions
Average Wear Time	<ul style="list-style-type: none"> • 24 hours – daily dressing required
Contra-Indications for Use	<ul style="list-style-type: none"> • Use with caution in the presence of significant hepatic or renal impairment. • Caution of use is required in patients known to be sensitive to systemic sulphonamides and in individuals known to have glucose-6-phosphate dehydrogenase deficiency. • Not recommended in heavily exudating wounds • Neonates
Application Guidance	A layer of Flamazine cream at least 3-5mm thick should be applied to the base of a previously cleansed wound or ulcer bed and cover with a suitable secondary dressing
PCOD Ordering information	<p>Order via pharmacy.</p> <p>The contents of one container are for the treatment of one person. 250g and 500g pots should be discarded 24 hours after opening. Tubes should be discarded 7 days after opening.</p>

Available Sizes	Via Pharmacy

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Podiatry Access	Foam
Product	BIATAIN
Company	Coloplast
Indications for Use	<ul style="list-style-type: none"> • Moderate to high exuding wounds • Leg ulcers, • Pressure ulcers • Non-infected diabetic foot ulcers. • May be used on patients who are in treatment for a local or systemic infection under the discretion of a physician
Average Wear Time	<ul style="list-style-type: none"> • Dependant on wound exudate. • May be left in place up to 7 days
Contra-Indications for Use	<ul style="list-style-type: none"> • Allergic reaction or known hypersensitivity to product
Application Guidance	Ensure that the plain side of foam is in contact (or facing the wound if using in conjunction with another dressing) with the wound and that the side with "Biatain" printed on it is facing away from the wound
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
5 x 7	141258
10 x 10	069682
10 x 20	3412
15 x 15	3413
20 x 20	3416

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	Silver Impregnated Foam
Product	BIATAIN AG
Company	Coloplast
Indications for Use	<ul style="list-style-type: none"> • Exuding leg ulcers, • Pressure sores, • Diabetic foot ulcers,
Average Wear Time	<ul style="list-style-type: none"> • Wear time can be up to 7 days. • depending on the amount of exudates
Contra-Indications for Use	<ul style="list-style-type: none"> • Previous reaction or known sensitivity to silver. • Not effective in dry wounds
Application Guidance	Ensure plain side of dressing is towards the wound and that the side printed with "Biatain Ag" is facing away from the wound. Non-adhesive may be secured with either tape or bandage.
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
Non Adhesive Range	
5 x 7	5105
10 x10	9622
10 x 20	8623
15 x 15	9625
20 x 20	9626

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	Foam
Product	ASKINA NON-ADHERENT FOAM including Heel shape
Company	B-Braun
Indications for Use	<ul style="list-style-type: none"> • Diabetic foot ulcers • Pressure ulcers • Leg ulcers
Average Wear Time	<ul style="list-style-type: none"> • Depends on exudate levels
Contra-Indications for Use	<ul style="list-style-type: none"> • Ulcers resulting from infections, such as tuberculosis, syphilis, deep fungal infections • Bites or third degree burns
Application Guidance	<p>Directly to wound as a primary dressing or as a secondary dressing. Remove backing and apply gently. Secure with bandage or Clinfast retention product.</p>
PCOD Ordering information	Product brand name, size and quantity required

Available Sizes	Procurement Codes
10 x10 cm	138128
10x20 cm	7241210
20 x 20	7242005
Heel 18.5 x 20.5 cm (225cm ²)	138111

Disclaimer: For more in depth product information please refer to company product leaflet.

TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION

Category Special Access	5-Layer Laminated Boot
Product	KERRABOOT
Company	Ark Therapeutics
Indications for Use	<ul style="list-style-type: none"> • Delayed-healing chronic wounds – leg ulcers, diabetic foot ulcers. • Can be used with necrotic areas in conjunction with vascular intervention. • If microvascular disease without PVD/PAD and necrosis, then it is safe to use as it will aid debridement to healthy tissue. • Can be used in patients with cellulitis and osteomyelitis in-conjunction with appropriate antibiotic therapy
Average Wear Time	<ul style="list-style-type: none"> • The product may increase levels of exudates initially which is normal. Kerraboot should be changed when the absorbent pad gels and becomes saturated with exudate. • This may be three times a day initially. Kerraboot would normally require to be changed daily in any case.
Contra-Indications for Use	<ul style="list-style-type: none"> • Kerraboot has not to be used along with any other dressings or compression bandages. • Discontinues use if rash evident. • Keep away from fire
Application Guidance	Easy to apply- slip boot over foot and leg and secure with Velcro strap just below the knee. This should be secure enough to prevent the boot slipping down the leg but not too tight. It can be adjusted throughout the day to accommodate any oedema.
PCOD Ordering information	Product brand name, size and quantity required

Shoe size & colour	Procurement Codes
7+ Extra Large Clear	3260239
< 3 Extra Small	3260221
7 Large	3013505
4-7 Small	3013497
Extra Large - White	3260247
Extra Small	3260270
Large	3223997
Small	3223989

Disclaimer: For more in depth product information please refer to company product leaflet.

**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

Company	Name	Contact Number	Products
Molnysky	Alan Gage	07796 675 620	Menilex Border
Advancis Medical	Lynsay Muir	07900 914 114	Activon
Smith & Nephew	Barbara Coventry	07968 726 522	Allevyn Gentle Range
			Allevyn
			Adoflex
			Adoflex
			Adoseb
Advanced Medical Solutions	Janey Mentipty	07791 742 957	Flamazine
			ActivHeal Hydrogel
			ActivHeal Hydrofiber
Systemics	Nora Haggerty	07895 806 093	Actifit
	Audrey Gillies	07771 373 485	Actisorb Cool
Aspen	Heather Guy	07720 293 769	Actisorb 220
			BioVita NA
			Promogram Prisma
B-Braun	Maureen Anderson		Askina Gel
			Askina Heel
3M	David Stevenson	07785318810	Tegaderm Foam Heel
		07795 240 121	Cavilon cream
Clinimed	Pam White	01628 0850100	Cavilon Film
			Cavilon Film
Coloplast	Catherine Patton	07801 233 767	Biatian Ag
Talley (Vacuum)	Jane Macdonald	07829 384 975	Biatian
			Venturi TNP
BioMonde	Katy Nicell	07896 687 010	Larvae/ BioFoam
	Ewan Murray	07989 258 758	
Crawford Healthcare Ltd	Duncan Griggs	07863 934 462	Flaminal Forte
			Flaminal Hydro
			Kerraboot
Convatec	Rosemary Reid	07801 857 082	Duoderm Extra Thin
			Aquacel
			Kaltostat
			Carboflex

Disclaimer: For more in depth product information please refer to company product leaflet.



**TISSUE VIABILITY SERVICES
WOUND MANAGEMENT FORMULARY
PRODUCT INFORMATION**

Uncontrolled when printed

Disclaimer: For more in depth product information please refer to company product leaflet.

Issue date May 2011
Review date May 2013

TVS